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Adapter Assembly on 12mm Electrode

Disassemble adapter and confirm your adapter fits a 12mm electrode.

Put adapter cap, spacer with chamfer pointer toward sensor end, O-rings (2), and lower body on to 12mm sensor. The o-rings will determine where the adapter is positioned on the sensor body.

Compress the adapter completely and screw the lower body to the upper adapter body.

Handling, Usage and Storage

When handling your pH probe, rinse the electrodes with distilled water before and after measuring a sample. Blot the end of the electrode with lint-free cloth to remove excess water. Never wipe the electrode to remove excess water – wiping can create static charges that interfere with the correct pH measurement.

Storage:

DO NOT STORE THE ELECTRODE IN DISTILLED OR DEIONIZED WATER. This will cause ions to leach out of the glass bulb and render the electrode useless. After storage, you may notice white KCl crystals deposited on your electrode. Such salt formation will not interfere with measurements. Simply rinse the electrode with distilled water to remove the crystals and blot dry before use.

When using a refillable pH electrode, the solutions should be filled up to, but not past, the refill hole. The refill hold must be open when measuring to ensure the fill solution flows properly through the reference junction.

pH electrodes are shipped with the pH bulb moist. Prior to using your electrode for the first time, follow these three steps to condition your electrode:

1. Remove the protective cap or boot from the bottom of the electrode and rinse the sensing glass with distilled or deionized water. (NOTE: Keep the protective cap of boot for use later during storage.)
2. Place the electrode in a beaker containing one of the liquids listed below (in order of ionic ability to condition the electrode). Soak for 20 minutes.
 - 4.0M KCl
 - 4.0 pH buffer
 - 7.0 pH buffer
3. After conditioning the electrode for 20 minutes, rinse the electrode with distilled or deionized water. The electrode is now ready for calibration and to measure pH.

Extend the Life of Your Electrode

When an electrode ages, it may exhibit sluggish or noisy readings. You can attempt to improve performance with the following procedures:

Reference problems:

A blocked reference junction is the most common problem of pH measurements. Symptoms include a slow response, off-scale readings and noisy readings. The procedure is specific to the type of electrode reference.

- GEL FILLED (non-refillable) – Soak the electrode in a beaker of warm water (60°C) for 15 minutes to remove the dried gel or salts from the junction. Then place in a beaker of warm 4M KCl solution. Set aside until it returns to room temperature. The gel should be moist, and the junction flow should be restored.

Glass Bulb Problems:

A glass membrane will get dirty over time. We suggest the following solutions for cleaning:

- PROTEIN – Use ASI Protein Cleaning Solution or wash in a solution of liquid soap, (about ½ teaspoon per 200mL warm water), using a soft cloth to gently wipe the pH glass – remember that pH glass is extremely delicate and breaks easily.
- INORGANIC SALTS – Wash in 0.1M HCl or EDTA (DO NOT SOAK); rinse with distilled water.
- GREASY FILMS – Wash in acetone or methanol (DO NOT SOAK); wash with liquid soap then rinse with distilled water.

After cleaning, place the electrode in a storage solution for 15 minutes prior to use. The best storage solution for a combination electrode is available from ASI in several convenient sizes. A pH 4 buffer is also acceptable.

pH Electrode Storage and Maintenance

To ensure a quick response and free-flowing liquid junction, the sensing element and reference junction must not be allowed to dry out. The following instructions apply to refillable electrodes. For gel-filled electrodes, consult instruction manual.

Routine Storage: Soak electrode in a pH Electrode Storage Solution. If a storage solution is unavailable, pH 4 buffer or 0.1 M KCl may be used.

Long-term Storage: Cover sensing element with its protective cap containing a several drops of pH/ORP storage solution. NOTE: Before returning to use, rinse well with deionized water, refill with appropriate fill solution and soak in pH 4 buffer for 1-2 hours.

Maintenance Cleaning your electrode between and after use will help extend the life of your electrode and avoid the cost of early replacement.

Routine Cleaning General: Soak electrode in 0.1 M HCl or HNO₃ for half an hour, followed by soaking it in storage solution for at least two hours. Reference junction - To clean a fouled junction: Heat electrode in 0.1 M KCl at 70°C for 10 minutes. Soak in storage solution for one hour prior to use.

Weekly Maintenance Inspect electrodes for scratches, cracks, salt crystal buildup, or membrane/junction deposits. Rinse off any salt buildup with distilled water and remove any membrane/junction deposits as directed in cleaning procedures below.


Removal of Membrane/Junction Deposits Grease and oil

- Rinse electrode with a mild soap Inorganic - Soak electrodes in 0.1M solution for 10 minutes Protein
- Soak in a 0.1M HCl solution for 10 minutes, rinse with DI water

After any of these cleaning procedures, soak the electrodes in storage solution overnight.

AUTO kit

The AUTO kit (pH measure and regulation) can be added to any EVO series unit. The following figures show the parts of the kit, plus an installation diagram.

Parts

- 1- EVO Equipment
- 2- Suction Tube (flexible)
- 3- Injection Tube (rigid)
- 4- Suction filter (place vertically on the bottom of the acid container)
- 5- pH Probe
- 6- Probe Holder
- 7- Flange **(not included)**
- 8- Injector (place with the arrow facing downwards)
- 9- Pipe nipple 3/8,1/2
- 10- Flange **(not included)**
- 11- Acid inlet (suction tube)
- 12- Acid outlet (injection tube)
- 13- pH probe connector (BNC)
- 14- pH4 calibration liquid
- 15- pH7 calibration liquid
- 16- Rubber cap for calibration